

# Lab Log

*Daniel Padfield*

**18/10/2017**

- Autoclaved the soil again
- Autoclaved 1 L of DI millicule water
- Autoclaved 1 L of 0.6% soft agar
- Autoclaved a spoon to weigh out soil
- Got some *SBW25φ2* from Floh. Around 1 mL at  $7.26 \times 10^8$ .
  - Grow up overnight
  - Transfer 60  $\mu$ L of bacteria and 10  $\mu$ L of phage (around  $7 \times 10^6$  pfus) into 6 mL of KB agar
  - Should give a concentration of around  $10^8$  phage/mL
  - Done this in triplicate
- Grow up *lacZ* and *WT* strains overnight. Should give concentration of around  $\sim 10^8$  cells in 60  $\mu$ L.
- Do these in triplicate
- Added 60  $\mu$ L of frozen overnight culture from first experiment (18/08/2017 *lacZ* and *WT*)

**Retrospectively work out density of the overnight stocks and phage**

**19/10/2017**

- Put 80g of soil into each 10cm x 10cm microcosm
  - Used autoclaved spoon
  - Placed scale in laminar flow hood (cleaned with ethanol before and after)
- Placed 5 mL ( $\sim 200$   $\mu$ L per microcosm) of *lacZ* and *WT* into separate 12 mL centrifuge tubes
  - Centrifuged for 15 minutes at max speed ( $\sim 4500$  r.p.m) on big centrifuge
  - Want to get to 5 mL per microcosm for inoculating ( $\sim 125$  mL in total)
  - Resuspended pellet into 2250  $\mu$ L, vortexed and placed 620  $\mu$ L, 620  $\mu$ L and 810  $\mu$ L into three different falcon tubes
  - Filled these three falcon tubes up to 40 mL, 40 mL and 45 mL respectively
    - \* This guaranteed the same concentration of sample in each falcon tube
  - Placed 5 mL of *lacZ* or *WT* strain into each microcosm
  - Froze ( $-80$  °C) 900  $\mu$ L of inoculate in 900  $\mu$ L of glycerol (25% final concentration)
- In the no phage treatments, we added 5 mL of M9
- Added 5 mL of phage to phage treatments
  - Place 900  $\mu$ L of bacteria + phage into three centrifuge tubes
  - Add 100  $\mu$ L (10%) chloroform into each tube (under fume hood)
  - Vortex rigorously
  - Centrifuge for 2-3 minutes at full speed (minifuge)
  - Take out supernatant and placed in a single tube (took out 800  $\mu$ L of each tube)
  - Put 40 mL of M9 into 6 tubes
  - Added 400  $\mu$ L into each tube ( $\sim 100$  fold dilution from the initial stock)
  - Shake each tube and add 5 mL into each microcosm
- Place microcosms into the 26 °C incubator (Level 1 incubator room)

**22/10/2017**

- setup 2 *WT* microcosms up for spot assays of phage
- used a crystal from  $T_0$
- autoclaved 2 L of KB agar. Put in small autoclave

**23/10/2017**

- phage spot test 10<sub>-1</sub> to 10<sub>-8</sub>
  - Floh's phage ( $7.26 \times 10^8$ )
  - Phage we inoculated with
  - Extra phage tube (should have lower concentration)
  - blank plate
- plated T<sub>0</sub> of WT and *lacZ P. fluorescens*

**24/10/2017**

- phage spot test was no good (streaky)
- likely that we did not wait long enough for the spot to dry
- set up more T<sub>0</sub> bacteria overnight

**25/10/2017**

- phage spot assay again
- plating only 4 spots per plate instead of 8
- counted T<sub>0</sub> counts (plated 30 µl)
  - *lacZ* @ 10<sup>-4</sup>: 18
  - WT @ 10<sup>-4</sup>: 64

**26/10/2017**

- check phage spot, 10<sup>-6</sup> looks to be the correct dilution
- further counts are done by putting the 10 µl phage into the 1% bacteria soft agar
  - vortex and plate
  - leave to dry
  - incubate overnight

**30/10/2017**

- grow up some T<sub>0</sub> bacteria overnight

**31/10/2017**

- serial dilution of phage to 10<sup>-6</sup>
- 50 µl bacteria + 10 µl phage + 5 ml soft agar
  - vortex
  - pour
- set samples up for inoculated phage, floh's phage stock and our phage dregs (each in triplicate)
- in 28 °C incubator overnight

**01/11/2017**

- sampled all soil microcosms (after 13 days)
  - ~2g of soil into a 12 mL centrifuge tube (with ~ 7 glass beads)
    - \* sample using the big end of a 5 mL pipette and weigh on scales under hood
  - add 10 mL of M9

- vortex for  $\sim 1$  minute
- add 900  $\mu\text{l}$  of sample to 900  $\mu\text{l}$  of 50% glycerol and freeze at  $-80^\circ\text{C}$  (bacteria + phage)
- make phage suspension for each tube
  - \* add 900  $\mu\text{l}$  of sample to centrifuge tube
  - \* add 100  $\mu\text{l}$  of chloroform under fume hood
  - \* vortex
  - \* centrifuge at full speed for 4 minutes
  - \* put supernatant into separate tube
    - samples 1-17 800  $\mu\text{l}$ , samples 18-48 750  $\mu\text{l}$
- names of samples
- 1-12: WT no phage
- 13-24: *lacZ* no phage
- 25-36: WT + phage
- 37-48: *lacZ* + phage
- Autoclaved things
  - 2 empty 500 mL bottles
  - M9 salts x10 (500 mL)
    - \* converted the weight of  $\text{Na}_2\text{HPO}_4$  to the amount of  $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$  needed.
    - \* MW of  $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$  is 268. MW of  $\text{Na}_2\text{HPO}_4$  is 142.
    - \* to convert from weight of  $\text{Na}_2\text{HPO}_4$  to  $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$  is to multiply by  $\frac{268}{142}$
  - 1 L of DI water
  - 2 L hard agar
  - centrifuge tubes.

sample	weight
1	2.1
2	2.0
3	1.9
4	2.0
5	2.0
6	2.0
7	2.0
8	2.0
9	2.1
10	1.9
11	2.2
12	2.1
13	2.1
14	1.9
15	2.2
16	2.2
17	1.9
18	2.2
19	2.0
20	2.3
21	2.3
22	2.2
23	1.9
24	2.1
25	2.1
26	2.2
27	1.9
28	2.0
29	2.2

sample	weight
30	2.3
31	2.0
32	2.2
33	2.2
34	2.2
35	2.0
36	2.0
37	2.3
38	1.9
39	2.2
40	1.9
41	2.2
42	2.0
43	2.0
44	2.2
45	2.1
46	2.2
47	2.2
48	2.1

**02/11/2017**

- set up phage spot tests of all soil phage suspensions against ancestral WT on soft agar.
  - just checking for presence/absence of phage
  - 5 mL per plate, 4 spots per plate, 12 spots in total
  - important to let the spot dry before moving
- poured plates (all Xgal)