

# Lab Log

Daniel Padfield

18/10/2017

- Autoclaved the soil again
- Autoclaved 1 L of DI millicule water
- Autoclaved 1 L of 0.6% soft agar
- Autoclaved a spoon to weigh out soil
- Got some *SBW25* from Floh. Around 1 mL at  $7.26 \times 10^8$ .
  - Grow up overnight
  - Transfer 60  $\mu$ L of bacteria and 10  $\mu$ L of phage (around  $7 \times 10^6$  pfus) into 6 mL of KB agar
  - Should give a concentration of around  $10^8$  phage/mL
  - Done this in triplicate
- Grow up *lacZ* and *WT* strains overnight. Should give concentration of around  $\sim 10^8$  cells in 60  $\mu$ L.
- Do these in triplicate
- Added 60  $\mu$ L of frozen overnight culture from first experiment (18/08/2017 *lacZ* and *WT*)

**Retrospectively work out density of the overnight stocks and phage**

19/10/2017

- Put 80g of soil into each 10cm x 10cm microcosm
  - Used autoclaved spoon
  - Placed scale in laminar flow hood (cleaned with ethanol before and after)
- Placed 5 mL ( $\sim 200$   $\mu$ L per microcosm) of *lacZ* and *WT* into separate 12 mL centrifuge tubes
  - Centrifuged for 15 minutes at max speed ( $\sim 4500$  r.p.m) on big centrifuge
  - Want to get to 5 mL per microcosm for inoculating ( $\sim 125$  mL in total)
  - Resuspended pellet into 2250  $\mu$ L, vortexed and placed 620  $\mu$ L, 620  $\mu$ L and 810  $\mu$ L into three different falcon tubes
  - Filled these three falcon tubes up to 40 mL, 40 mL and 45 mL respectively
    - \* This guaranteed the same concentration of sample in each falcon tube
  - Placed 5 mL of *lacZ* or *WT* strain into each microcosm
  - Froze ( $-80$   $^{\circ}$ C) 900  $\mu$ L of inoculate in 900  $\mu$ L of glycerol (25% final concentration)
- In the no phage treatments, we added 5 mL of M9
- Added 5 mL of phage to phage treatments
  - Place 900  $\mu$ L of bacteria + phage into three centrifuge tubes
  - Add 100  $\mu$ L (10%) chloroform into each tube (under fume hood)
  - Vortex rigorously
  - Centrifuge for 2-3 minutes at full speed (minifuge)
  - Take out supernatant and placed in a single tube (took out 800  $\mu$ L of each tube)
  - Put 40 mL of M9 into 6 tubes
  - Added 400  $\mu$ L into each tube ( $\sim 100$  fold dilution from the initial stock)
  - Shake each tube and add 5 mL into each microcosm
- Place microcosms into the 26  $^{\circ}$ C incubator (Level 1 incubator room)

22/10/2017

- setup 2 *WT* microcosms up for spot assays of phage
- used a crystal from  $T_0$
- autoclaved 2 L of KB agar. Put in small autoclave

**23/10/2017**

- phage spot test 10<sub>-1</sub> to 10<sub>-8</sub>
  - Floh's phage ( $7.26 \times 10^8$ )
  - Phage we inoculated with
  - Extra phage tube (should have lower concentration)
  - blank plate
- plated T<sub>0</sub> of WT and *lacZ P. fluorescens*

**24/10/2017**

- phage spot test was no good (streaky)
- likely that we did not wait long enough for the spot to dry
- set up more T<sub>0</sub> bacteria overnight

**25/10/2017**

- phage spot assay again
- plating only 4 spots per plate instead of 8
- counted T<sub>0</sub> counts (plated 30 µl)
  - *lacZ* @ 10<sup>-4</sup>: 18
  - WT @ 10<sup>-4</sup>: 64

**26/10/2017**

- check phage spot, 10<sup>-6</sup> looks to be the correct dilution
- further counts are done by putting the 10 µl phage into the 1% bacteria soft agar
  - vortex and plate
  - leave to dry
  - incubate overnight

**30/10/2017**

- grow up some T<sub>0</sub> bacteria overnight

**31/10/2017**

- serial dilution of phage to 10<sup>-6</sup>
- 50 µl bacteria + 10 µl phage + 5 ml soft agar
  - vortex
  - pour
- set samples up for inoculated phage, floh's phage stock and our phage dregs (each in triplicate)
- in 28 °C incubator overnight

**01/11/2017**

- looked at phage spots (10<sup>-6</sup> dilution)
- calculated how many phage are there

$$\text{Phage concentration} = \frac{\text{number of plaques}}{\text{dilution factor} \times \text{volume added}}$$

Volume added is almost always 10 µl. I wrote a mini function to calculate this

```
d <- data.frame(dregs = c(3,4,6), phage = c(6,2,4))

# back calculate number in there
# PFU/ml = plaque number / (dil fac * volume added)
plaq_num <- function(x, dilfac){return(x/(dilfac*0.01))}

d <- dplyr::mutate_at(d, c('dregs', 'phage'), plaq_num, dilfac = 10^-6)

knitr::kable(d)
```

dregs	phage
3e+08	6e+08
4e+08	2e+08
6e+08	4e+08

```
# concentration of phage added
mean(d$phage)
```

```
## [1] 4e+08
```

```
# concentration of other phage
mean(d$dregs)
```

```
## [1] 433333333
```

- sampled all soil microcosms (after 13 days)
  - ~2g of soil into a 12 mL centrifuge tube (with ~ 7 glass beads)
    - \* sample using the big end of a 5 mL pipette and weigh on scales under hood
  - add 10 mL of M9
  - vortex for ~ 1 minute
  - add 900 µl of sample to 900 µl of 50% glycerol and freeze at -80 °C (bacteria + phage)
  - make phage suspension for each tube
    - \* add 900 µl of sample to centrifuge tube
    - \* add 100 µl of chloroform under fume hood
    - \* vortex
    - \* centrifuge at full speed for 4 minutes
    - \* put supernatant into separate tube
      - samples 1-17 800 µl, samples 18-48 750 µl
- names of samples
- 1-12: WT no phage
- 13-24: *lacZ* no phage
- 25-36: WT + phage
- 37-48: *lacZ* + phage
- Autoclaved things
  - 2 empty 500 mL bottles
  - M9 salts x10 (500 mL)
    - \* converted the weight of  $Na_2HPO_4$  to the amount of  $Na_2HPO_4 \cdot 7H_2O$  needed.
    - \* MW of  $Na_2HPO_4 \cdot 7H_2O$  is 268. MW of  $Na_2HPO_4$  is 142.
    - \* to convert from weight of  $Na_2HPO_4$  to  $Na_2HPO_4 \cdot 7H_2O$  is to multiply by  $\frac{268}{142}$

- 1 L of DI water
- 2 L hard agar
- centrifuge tubes.

sample	weight_g
1	2.1
2	2.0
3	1.9
4	2.0
5	2.0
6	2.0
7	2.0
8	2.0
9	2.1
10	1.9
11	2.2
12	2.1
13	2.1
14	1.9
15	2.2
16	2.2
17	1.9
18	2.2
19	2.0
20	2.3
21	2.3
22	2.2
23	1.9
24	2.1
25	2.1
26	2.2
27	1.9
28	2.0
29	2.2
30	2.3
31	2.0
32	2.2
33	2.2
34	2.2
35	2.0
36	2.0
37	2.3
38	1.9
39	2.2
40	1.9
41	2.2
42	2.0
43	2.0
44	2.2
45	2.1
46	2.2
47	2.2
48	2.1

**02/11/2017**

- set up phage spot tests of all soil phage suspensions against ancestral WT on soft agar.
  - just checking for presence/absence of phage
  - 5 mL per plate, 4 spots per plate, 12 spots in total
  - important to let the spot dry before moving
- poured plates (all Xgal)